



Simple, Economical Solution for Lighter Weight Chairs and Scooters

Affordable lift transports and stores a power chair or scooter inside a cargo area. Single axis lift rotates manually with powered lifting; arm height adjusts to application.

Features & Benefits

- Adjustable arm height
- Manual rotation
- Weather-resistant hand control
- Passenger side standard; driver side available
- Offset post allows the mobility device to be safely unloaded on the side of the vehicle
- Docking device included based on application

| Specifications | AL215 Axis I Lift |
|------------------------|-------------------------|
| Accommodates | Power chairs & scooters |
| Lifting Capacity | 250 lbs. (113 kg) |
| Installed Weight | 42 lbs. (19 kg) |
| Adjustable Boom Height | 30" – 40" |
| Adjustable Boom Length | 21.25" – 32.35" |
| Adjustable Boom Pitch | 5° |
| Securement | Docking device included |
| Options | Battery pack |
| Warranty | 3-year transferable |

Key Features



Easy-to-use hand control with Soft-Touch buttons.



Extensive Docking Device selection (DD19 & DD17 shown)