


Take Yourself Higher

Residential Vertical Platform Lift



 **BRUNO**[®]
for your independence



MADE IN AMERICA

Learn more at www.bruno.com



RESIDENTIAL VPL

Count on Bruno Quality

- **EASY-TO-USE** — Pressure control rocker switch up/down operation. Automatic self-lowering ramp folds down for smooth roll on/off access; folds up when in use to create a safety barrier.
- **HIGH-PERFORMANCE** — Lift capacity of 750 lb/340 kg. Smooth, quiet ride. Made-in-the-USA quality.
- **SAFE** — Unit shuts down when obstruction is touched under platform. Emergency stop switch on platform.
- **RELIABLE** — Continuously charged DC battery-powered motor on lifts 10-feet and above ensures operation even in a power outage. Affordable AC in 4- and 6-foot units (DC battery powered optional) Gold warranty: 2 years, major components; 1 year, parts

Fit Your Space and Life

- **Minimum Space. Maximum Reach** — Space-efficient Bruno wheelchair lifts can access up to 14 feet.
- **Indoor or Outdoor** — Weather-resistant controls and specialized coating protect the unit. An optional canopy can protect you as well. Combat extreme winters with the optional cold weather package.



Options to Meet Your Needs

- Top landing gate
- Landing call/send station
- Platform gate
- Platform canopy
- Multiple platform sizes and configurations
- Cold-weather package
- Paddle controls
- DC battery powered in 4- and 6-foot units



Give Your Life a Lift

Gain access to your porch, deck, or floor-to-floor without sacrificing space with a compact Bruno residential vertical platform lift. The quiet, smooth lift accommodates any wheelchair, power chair or scooter. Give yourself the freedom to get in and out of your home with just a push of a button. You're in control.

Learn more at www.bruno.com





Innovative solutions for your indoor and outdoor lifestyle



Indoor/Outdoor Stairlifts



Vertical Platform Lifts



Scooter/Wheelchair Lifts



Valet® Signature Seating

About Bruno

Bruno Independent Living Aids has helped improve the lives of people with limited mobility for over 30 years. Engineered and manufactured in the USA, Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/powerchair lifts help people easily transport their mobility devices. In addition, Bruno's Valet Signature Seating is the world leader in revolutionary turning automotive seats to assist people in getting in and out of their vehicles more easily. A family-owned business, Bruno takes pride in their unwavering focus on quality and safety and their position as industry innovation leader.

Proudly made in the U.S.A.



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2016

MKT-VPL-TT