



STAY IN THE HOME
YOU LOVE

www.bruno.com



BRUNO'S BEST-SELLING STAIRLIFT, THE ELAN

Bruno's most popular stairlift can help restore full access to your home. Designed for straight indoor staircases, the Elan offers a mix of high-quality standard functionality, unsurpassed reliability, and ease-of-use design.

With Made-in-the-USA dependability, Bruno helps you stay in your home with comfort and ease.

Learn more at www.bruno.com.



101 Mobility®

Experts in Mobility and Accessibility Solutions

101 Mobility has people who are passionate in providing the optimal solutions and service for their customers. From the first phone call to the final installation, we pride ourselves on unconditional support – no matter what.

www.101mobility.com

- Locally owned and operated
- Outstanding service and support
- Expert installation
- Extended 101 warranty protection



ELAN FEATURES


- Armrest rocker switch for easy operation
- Anodized rail with covered gear rack gives a sleek appearance
- Vertical rail design enables tight installation to within 5 in (127 mm) of wall
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Continuous charge strips along rail recharge unit wherever it stops
- Two wireless call/sends allow unit to be controlled remotely
- Two 12v batteries ensure dependable performance even in a power outage
- Soft start/stop gives a smooth ride from start to finish
- Obstruction safety sensors, seat belt ensure safety
- Onboard audio diagnostics for easy service




Padded, generous seat.



ELAN OPTIONS



Power or manual folding rail for narrow hallway or when doorway is at the bottom of stairs. Manual operation or push button automatic.



Power swivel seat for effortless exit. Armrest switch control or wireless call/send.



Power folding footrest automatically flips up/down when seat is raised/lowered. Arm switch control optional.



Offset swivel seat makes getting on/off easy.



Sensors ensure safety.



Two wireless call/sends allow remote control access.



Compact, flip-up design for plenty of open space on stairway.

“Staff was very professional and helpful. The chair is wonderful and I will recommend the product and company.”

V. KOPS - CA



“Works smoothly and quietly to move my Mother from floor to floor. Appreciate the safety features as well as the ease of operation”

C. Savan - OH



TECHNICAL SPECIFICATIONS

Seat Pad Size: 17.5 in wide x 15 in deep
(445 mm wide x 381 mm deep)

Width Between Arms: 20 in (508 mm)

Footrest Size: 10.5 in long x 11 in wide
(267 mm long x 279 mm wide)

Vertical Rail:

Installs within 5 in (127 mm) of the wall

About 101 Mobility

101 Mobility is your local expert in mobility and accessibility solutions. Independently owned and operated, each 101 Mobility location sells, installs and services stairlifts, platform lifts, wheelchair ramps and patient lifts to adapt to you and your family's changing needs.

About Bruno

Family-owned Bruno Independent Living Aids has improved the lives of people with limited mobility for over 30 years. Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/ powerchair lifts help people easily transport their mobility devices. In addition, Bruno's Valet Signature Seating is the world leader in revolutionary turning automotive seats to assist people in getting in and out of their vehicles more easily.

Proudly made in the USA.



ACCREDITED BUSINESS



TRUSTPILOT | CONSUMER AFFAIRS

101 Mobility[®]
Experts in Mobility and Accessibility Solutions

 **BRUNO**[®]