



### Folding Universal Sliding Transfer Bench

### Bariatric Transfer Bench Tool-free Removable Back and Arm

Item #	UOM
RTL12075	2/cs

Item #	Product Description	UOM
12025KD-2	Knocked Down, Tool-free Removable Back and Arm	2/cs
12025KD-1	Knocked Down, Tool-free Removable Back and Arm	1/cs

- Transfer Bench slides left and right along frame for easy entry into and out of tub.
- Seat height adjusts in 1/2" increments.
- Folds flat for transport and storage. (Figure A)
- Comes standard with removable soap dish. (Figure B)
- Comes with safety net to catch dropped personal items and prevent bending.
- Comes with seat belt for added safety.
- Retail packaged. (Figure D)

- 500 lb. weight capacity.
- Tool-free assembly of back, legs and arm.
- A-frame provides additional stability.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- All transfer benches are reversible to accommodate any bathroom.
- Euro-style clip ensures a secure fit. (Figure A)
- No exposed hardware to injure patient.
- Large, locking suction cups for added safety. (Figure B)

specifications	
Seat Dimensions	22.5" (W) x 18.5" (D)
Height	21" - 25"
Weight	16 lbs.
Carton Shipping Weight	38 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

specifications	
Outside Legs	32.5" (W) x 24" (D)
Seat Dimensions	26.5" (W) x 16.5" (D)
Seat Height	18.5" - 23.5"
Weight	14 lbs.
Carton Shipping Weight (-1   -2)	14 lbs.   27 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime