

# AL160

## PROFILE SCOOTER LIFT



### SEE THE VIEW, NOT THE LIFT

The open platform design of the AL160 creates a simple, streamlined silhouette. Its understated presence on a vehicle is stylish and ideal for compact cars.

Harmar's AL160 Profile Scooter Lift safely and securely transports power scooters with a wheelbase between 30" and 42" – no scooter modifications needed. The automatic hold-down foot secures your scooter automatically with the push of a single button. Even larger power scooters can be transported with this profile lift – up to 350 lb.

Single switch operation and key lock



Hold-down foot



*The AL160 profile lift is ideal for my scooter. Its open platform frame allows for great visibility when driving without my scooter and has plenty of capacity for when I need to take it along.*

Rear wheel cradle

## AL160 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	350 lb	<b>ACCOMMODATES</b>	Scoters up to 42" wheelbase
<b>INSTALLED WEIGHT</b>	74 lb	<b>PLATFORM SIZE</b>	25" x 25" adjustable cradle extends up to 33"
<b>HITCH CLASS</b>	I or greater	<b>POWERED</b>	Fully powered
<b>WARRANTY</b>	3-year	<b>INCLUDED</b>	AL123 class III hitch adaptor, AL107 hitch post
<b>AVAILABLE OPTION(S)</b>	<ul style="list-style-type: none"> <li>- Battery pack</li> <li>- Swing-away</li> <li>- Scooter cover</li> </ul>	<b>SECUREMENT</b>	<ul style="list-style-type: none"> <li>- Hold-down foot</li> <li>- Rear wheel cradle</li> </ul>