

Side Door Hybrid Van Lift AL690



The Industry's Only Side Load Hybrid

Convenient for drive on and off side loading. The AL690 lifts and loads both wheelchairs and scooters securely into vans and most mini-vans.

Features & Benefits

- Compatible with many power chairs and shorter scooters
- Drive on/Drive off platform
- No scooter or chair modifications required
- No drilling required on most applications
- Connects to existing second row seat hardware
- Third row seating spared in most applications
- Remote hand control (wireless option available)
- 275 lbs. (125 kg) capacity



HARMAR EXCLUSIVE



Specifications	AL690 Side Door Hybrid Platform Lift
Accommodates	Power Chairs & Scooters
Lifting Capacity	275 lbs. (125 kg)
Installed Weight	186 lbs. (84 kg)
Platform Size	17" x 28" adjustable width; 38" length
Securement	Manual Straps
Options	Self-tensioning straps, battery pack, wireless remote
Warranty	3-year transferable



800-833-0478 | www.harmar.com

