



SECURITY POLE™




FEATURES



ERGONOMIC GRIP



SUPPORTS UP TO 300 LBS




TENSION MOUNTED - NO SCREWS OR DRILLING



ADJUSTS TO FIT 7' - 10' CEILINGS



USE IN ANY ROOM



MAKES SITTING AND STANDING EASY

SPECIFICATIONS

- *Materials: Rust resistant zinc plated steel with durable powder coat finish*
- *Fits Ceiling Heights: Range of 7'-10' ft.*
- *Weight of Product: 20 lbs*
- *Package Dimensions: 54.00" L x 12.50" W x 3.25" H; 22 lbs*
- *Case of 2 Dimensions: 54.00" L x 12.50" W x 6.5" H; 45 lbs*

COMES IN TWO COLORS

 BLACK #1150-B

 WHITE #1150-W

PRODUCT #1150