

Options for Liko Viking M, Liko Viking L and Liko Viking XL:

- **Armrest** – Retractable armrests attach onto lift for patient gait training and user maneuverability of the lift
- **Bag for SlingBars** – Attaches on mobile lift to provide storage for SlingBars
- **Battery** – Lift battery for spare usage
- **Cable reel EU** – Attaches on mobile lift to enable mounting of external charge and charger cable
- **LikoScale 350** – Scale can weigh patients up to 772 lbs (350 kg)
- **Quick-Release Hook** – Attaches to Liko sling bars not already equipped with Quick-Release Hook for safe and easy changing of accessories
- **SlingBar Cover Paddy 30** – Covers sling bar to create padding for patients who risk hitting their heads on sling bar due to involuntary movement. Compatible with:
 - *Universal SlingBar 350*
 - *Universal SlingBar 450*
 - *Universal SlingBar 600*
- **Slings** – A variety of Liko's slings are compatible with Multirall and its sling bars to achieve safe, versatile and comfortable patient suspension and transfer. With over 20 different sling models, sized from XS to bariatric, and many different materials, Liko slings can be used for:
 - General lifting situations
 - Hygiene lifting in connection to toiling or bathing
 - Repositioning
 - Gait training and rehabilitation
 - Bariatric lifting
 - Specialized lifting situations
- **Sling bars** – Eight recommended sling bars for versatile and comfortable patient suspension:
 - **SlingBar Mini 220** – Weight capacity: 450 lbs (205 kg)
 - **Sling Cross-bar 450** – Weight capacity: 660 lbs (300 kg)
 - **Sling Cross-bar 670** – Weight capacity: 660 lbs (300 kg)
 - **Universal SideBars 450 including bag** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 350** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 450** – Weight capacity: 660 lbs (300 kg)
 - **Universal SlingBar 600** – Weight capacity 660 lbs (300 kg)
 - **Universal TwinBar 670 Twin** – Weight capacity: 660 lbs (300 kg)
- **Stretchers** – One recommended stretchers to enable easy horizontal lifting:
 - **OctoStretch Lift Stretcher**
 - Weight capacity: 440 lbs (200 kg)

- Unit dimensions:
 - 63" x 20" (1600 mm x 500 mm)
 - 67" x 35.4" (1700 mm x 900 mm)
- **Wall-mounted Battery Charger** – Mounts to the wall for easy battery charging

Options for Liko Viking L and Liko Viking XL only:

- **Stretchers** – Three recommended stretchers to enable easy horizontal lifting:
 - **FlexoStretch**
 - Weight capacity: 660 lbs (300 kg)
 - Unit width (adjustable): 27.6" – 34.4" (700 mm – 950 mm)
 - Unit length:
 - 65" (1650 mm)
 - 73" (1850 mm)
 - **LikoStretch Mod 600 IC Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit width:
 - 28" (711 mm)
 - 32" (813 mm)
 - **LikoStretch 1900 Lift Stretcher**
 - Weight capacity: 550 lbs (250 kg)
 - Unit dimensions: 74.80" x 35.79" x 23.35" (1900 mm x 909 mm x 593 mm)
- **Wheel Set** – Changes height of the lift's base to 4" (102 mm) for optimum customization